



**GeteinAnimal  
Medical**



**A Compact and Automated  
PCR Lab for Pet Pathogens**

**GN 7000 Vet**  
Akso PCR System

Getein Animal Medical Technology Co., Ltd.

# GN 7000 Vet

## Akso PCR System



### 1 Sample In - Result Out

- Core innovative technology
- Turnaround time as fast as  $\leq 35$  minutes

### 2 High Throughput and Quick Test

- Dual module design and 4 channels for simultaneous testing of multiple items.
- Special heated base with microfluidic reagent kit can realize fast heating and cooling.

### 3 Compact, Portable and Easy to Operate

- No need for multiple sets of equipment and cumbersome operational processes.

### 4 Low Cost and Easy Storage

- Molded consumables
- Simple assembly process
- The lyophilized reagent is embedded inside the kit and can be transported at room temperature.




More Than Just Quick

## Technical Parameter

Screen size: 10.1 inch  
Dimensions: 419 mm(L) x 310 (W) x 233 mm (H)  
Database storage volume: 100,000 pieces  
Working environment: 10-30 °C  
Weight: 9.5 kg



Professional testing report



# Test Items



:12 Pathogens ★:Test Items Available

Test Mode	Pathogen	Specimen
Respiratory Diseases		
Single Test	★CDV	Mouth, Nose and Conjunctiva Swab
Single Test	CPIV	Mouth, Nose and Conjunctiva Swab
Single Test	CAV- I	Mouth, Nose and Conjunctiva Swab
Single Test	CAV- II	Mouth, Nose and Conjunctiva Swab
Single Test	Bb	Mouth, Nose and Conjunctiva Swab
Combo Test	★ CDV/CPIV/CAV/Bb	Mouth, Nose and Conjunctiva Swab
Digestive Tract Diseases		
Single Test	CCV	Feces or Anal Swab
Single Test	★ CPV	Feces or Anal Swab
Combo Test	★ CCV/CPV	Feces or Anal Swab
Routine Physical Examination		
Combo Test	★ CDV/CCV/CPV/CAV	Mouth, Nose and Conjunctiva Swab+Feces or Anal Swab
Anemia		
Single Test	Bab	EDTA Anticoagulated Whole Blood
Single Test	Bgi	EDTA Anticoagulated Whole Blood
Single Test	APH	EDTA Anticoagulated Whole Blood
Single Test	Ehr	EDTA Anticoagulated Whole Blood
Combo Test	★ Bab/Bgi	EDTA Anticoagulated Whole Blood
Combo Test	APH/Ehr	EDTA Anticoagulated Whole Blood
Combo Test	★ Bab/Bgi/APH/Ehr	EDTA Anticoagulated Whole Blood
Zoonotic Diseases		
Single Test	★ TOX	EDTA Anticoagulated Whole Blood
Single Test	Lep	EDTA Anticoagulated Whole Blood
Combo Test	★ TOX/Lep	EDTA Anticoagulated Whole Blood



: 13 Pathogens ★ :Test Items Available

Test Mode	Pathogen	Specimen
Respiratory Diseases		
Single Test	FHV	Mouth , Nose and Conjunctiva Swab
Single Test	FCV	Mouth , Nose and Conjunctiva Swab
Single Test	MF	Mouth , Nose and Conjunctiva Swab
Single Test	CF	Mouth , Nose and Conjunctiva Swab
Single Test	Bb	Mouth , Nose and Conjunctiva Swab
Combo Test	★ FHV/FCV/MF/CF/Bb	Mouth , Nose and Conjunctiva Swab
Digestive Tract Diseases		
Single Test	FCoV	Feces or Anal Swab
Single Test	★ FPV	Feces or Anal Swab
Single Test	TF	Feces or Anal Swab
Single Test	GL	Feces or Anal Swab
Combo Test	★ FCoV/FPV/TF/GL	Feces or Anal Swab
Feline Anemia		
Single Test	FeLV	EDTA Anticoagulated Whole Blood
Single Test	FIV	EDTA Anticoagulated Whole Blood
Combo Test	★ FeLV/FIV	EDTA Anticoagulated Whole Blood
Zoonotic Diseases		
Single Test	★ TOX	EDTA Anticoagulated whole Blood
Single Test	Bh	EDTA Anticoagulated Whole Blood
Combo Test	★ TOX/Bh	EDTA Anticoagulated Whole Blood
Feline Stomatitis		
Combo Test	★ FHV/FCV/FeLV/FIV	Mouth, Nose and Conjunctiva Swab
Routine Physical Examination		
Combo Test	★ FPV/FCoV/FHV/FCV	Mouth, Nose and Conjunctiva Swab+Feces or Anal Swab



**GeteinAnimal Medical**

Add.: No. 9 Bofu Road, Luhe District, Nanjing, 211505, China

Tel.: +86-25-68568508/68568594 Fax: +86-25-68568500

E-mail: sales@getein.com.cn; overseas@getein.com.cn

Web.: www.getein.com

